


REKODSAYA PORTAL - PREREQUISITE CHECK IN

Table of contents

1. VC Appointment Blink Button	3
2. Appointment (Check-In)	4
3. Consent Checklist	5
4. Patient Instructions	6
5. Self-monitoring - Vital Sign	7
6. Self-monitoring - Final Vital Sign	8
7. Self-monitoring - BMI	9
8. Self-monitoring - Final BMI	10
9. Self-monitoring - Glucose	11
10. Self-monitoring - Final Glucose	12
11. Submit Prerequisites	13
12. Submit Prerequisites (Continued)	14
13. Info (Self-monitoring)	15
14. Self-monitoring - Observation	16
15. Self-monitoring - Observation (Select Panel)	17
16. Self-monitoring - Observation (Vital Signs)	18
17. Self-monitoring - Observation (Glucose)	19
18. Self-monitoring - Observation (Anthropometry)	20

1. VC Appointment Blink Button

1. System displays the Rekod Saya Homepage.
2. Click on the [VC APPOINTMENT] button to open appointment page.



2. Appointment (Check-In)

1. Click on the check-in button to fill in the pre-requisite form.

Note:

The button will display one (1) hour before actual appointment time.

The screenshot displays the REKODSAYA PORTAL interface. The top navigation bar includes the Malaysian coat of arms, home, health record, blood donations, and health centre links. A 'VC APPOINTMENT' button is visible in the top right. The main content area is titled 'Appointment' and 'Booking'. It features three tables: 'BOOKED', 'FULFILLED', and 'CANCELLED/EXPIRED'. The 'BOOKED' table has one entry with a 'CHECK-IN' button and a '1' hour timer. The 'FULFILLED' table is empty. The 'CANCELLED/EXPIRED' table has one entry with a 'CANCELLED' status.

BOOKED									
No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	Option
1	HTAN - KLINIK VC INDIVIDUAL	25-Aug-2023	16:00 - 17:00	BILIK SAKIT 1	Addiction & Substance Abuse	Virtual	Hospital Tuanku Ampuan Najihah, Kuala Pilah	BOOKED	CHECK-IN 1 CANCEL

FULFILLED									
No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	
No data present in this section.									

CANCELLED/EXPIRED									
No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	
1	HTAN - KLINIK VC INDIVIDUAL	25-Aug-2023	16:00 - 17:00	BILIK SAKIT 1	Addiction & Substance Abuse	Virtual	Hospital Tuanku Ampuan Najihah, Kuala Pilah	CANCELLED	

RekodSaya Portal © 2023 Session Timeout in 14 Min 56 Sec

3. Consent Checklist

1. Click on [Next] button to proceed next page.

Note:

Read the consent before proceed to next page.

The screenshot displays the REKODSAYA PORTAL interface. At the top, there is a navigation bar with the Malaysian coat of arms, a home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. On the right side of the top bar, there is a 'VC APPOINTMENT' button, a language selector set to 'EN', a notification icon with a '3' badge, and the name 'KASIH BALQIS'. Below the navigation bar is a 'MENU' sidebar with various health-related icons and labels: Health Events, Appointment (highlighted), Allergy, Imaging, Family History, Comorbidities, Medication, Laboratory, Observation, Immunization, Transfusion History, Referral, Self Assessment, and Billing. The main content area is titled 'Appointment' and 'Booking'. A progress indicator shows four steps: 1. Consent (active), 2. Patient Instructions, 3. Self-monitoring, and 4. Done. Below the progress indicator is a 'Consent Checklist' section containing a text box with the following text: 'Virtual Consultation Service is the delivery of virtual (virtual), direct (live) and interactive health services that include clinical consultation and customer treatment plans between health personnel and customers. This service is an initiative implemented to ensure social distancing and provide continuous health services at par with technological advances. I have been explained by the Medical Practitioner about the implementation of this Virtual Consultation Service and agree;'. Below this text is a numbered list of seven terms and conditions. At the bottom of the consent box, there are four buttons: 'Cancel', 'Back', 'Next', and a circular button with the number '1'. At the very bottom of the page, there is a footer with 'RekodSaya Portal © 2023' and a 'Session Timeout in 14 Min 59 Sec' warning.

REKODSAYA PORTAL - PREREQUISITE CHECK IN

4. Patient Instructions

1. Key in Comment field.
2. Click on [Next] button to proceed next page.

The screenshot displays the REKODSAYA PORTAL interface. At the top, there is a navigation bar with the logo, home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. A red 'VC APPOINTMENT' button is visible on the right. Below the navigation bar, a 'MENU' sidebar lists various health services: Health Events, Appointment (highlighted), Allergy, Imaging, Family History, Comorbidities, Medication, Laboratory, Observation, Immunization, Transfusion History, Referral, Self Assessment, and Billing. The main content area is titled 'Appointment' and 'Booking'. A warning message reads: 'If you feel unwell, kindly see your health care provider physically for a more thorough exam'. The 'PREREQUISITES FORM' is shown with a progress bar indicating four steps: 1. Consent, 2. Patient Instructions (current step), 3. Self-monitoring, and 4. Done. The 'Patient Instructions' section contains an 'Instructions' field with the text 'eat' and a 'Comment' field with the text 'ok'. A blue circle with the number '1' is positioned over the 'Comment' field. Below the form are 'Cancel' and 'Next' buttons, with a blue circle containing the number '2' over the 'Next' button. At the bottom left, it says 'RekodSaya Portal © 2023' and at the bottom center, a purple bar indicates 'Session Timeout In 14 Min 58 Sec'.

5. Self-monitoring - Vital Sign


1. Key in Vital Sign field.
2. Click on [SAVE] button.

The screenshot displays the REKODSAYA PORTAL - PREREQUISITE CHECK IN interface. The top navigation bar includes the portal logo, a home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. A red 'VC APPOINTMENT' button is visible on the right, along with language and user settings icons. The left sidebar contains a 'MENU' with various health-related categories, with 'Appointment' highlighted. The main content area is titled 'Appointment' and features a 'VITAL SIGN PANEL' with a blue circle containing the number '1'. The panel contains five input fields with their respective values and normal ranges: Systolic Blood Pressure (66 mmHg, Normal Range: 120 - 90 mmHg), Diastolic Blood Pressure (55 mmHg, Normal Range: 80 - 60 mmHg), Heart Rate (99 /min, Normal Range: 60 - 100 beats per minute), Body Temperature (37°C, Normal Range: 35°C - 37.5°C), and Oxygen Saturation (89%, Normal Range: 95% - 100%). The Respiratory Rate (67 /min, Normal Range: 12 - 20 breaths per minute) field is highlighted with a green border. At the bottom of the panel, there is a red 'SAVE' button with a white plus icon and a blue circle containing the number '2'. The footer of the page shows 'RekodSaya Portal © 2023' and a purple 'Session Timeout in 14 Min 57 Sec' warning.

REKODSAYA PORTAL - PREREQUISITE CHECK IN

6. Self-monitoring - Final Vital Sign

1. The system displays pop-up message.
2. Click on [OK] button.



7. Self-monitoring - BMI

1. Key in BMI questions field.
2. Click on [SAVE] button.

The screenshot displays the 'Self-monitoring' section of the REKODSAYA PORTAL. The interface includes a top navigation bar with 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. A red 'VC APPOINTMENT' button is located in the top right. The main content area is divided into 'Appointment' and 'Booking' tabs. The 'Appointment' tab is active, showing a progress bar with four steps: 'Consent', 'Patient Instructions', 'Self-monitoring' (current step), and 'Done'. Below the progress bar, a message states: 'Self-monitoring is not mandatory to proceed with the virtual consultation session'. The 'Vital Sign' section is expanded to show the 'BMI' panel. The BMI panel contains the following fields and values:

Field	Value
Height (cm)	150
Weight (kg)	55
BMI	24.44
Waist Circumference (cm)	42

The BMI panel also displays a reference range: 'Normal (BMI 20 - 24)'. A red 'SAVE' button is located at the bottom of the BMI panel, highlighted with a red circle and the number 2. A blue circle with the number 1 is placed over the 'BMI PANEL' header. The 'GLUCOSE' section is partially visible below the BMI panel. The footer of the page shows 'RekodSaya Portal © 2023' and a 'Session Timeout in 14 Min 58 Sec' warning.

8. Self-monitoring - Final BMI

1. The system displays pop-up message.
2. Click on [OK] button.

The screenshot displays the RekodSaya Portal interface. The top navigation bar includes the national emblem, 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. A red 'VC APPOINTMENT' button is visible on the right. The left sidebar contains a 'MENU' with various health-related options. The main content area is titled 'Appointment' and 'Booking', with progress indicators for 'Consent', 'Patient Instructions', 'Self-monitoring', and 'Done'. A 'Self-monitoring is not mandatory to proceed with the virtual consultation session' message is shown. The 'Vital Sign' section includes a 'BMI' panel with input fields for Height (cm) (150), Weight (kg) (55), and Waist Circumference (cm) (42). The calculated BMI is 24.44, with a normal range of 20-24. A 'SAVE' button is present. A pop-up dialog box is overlaid on the BMI panel, containing a warning icon, the text 'FINAL BMI 1', and the question 'Are you sure you want to SAVE this BMI? This action is not be able to update once saved'. The dialog has 'Cancel' and 'OK' buttons, with a blue circle '2' next to the 'OK' button.

9. Self-monitoring - Glucose

1. Key in Glucose questions field.
2. Click on [SAVE] button.

The screenshot displays the REKODSAYA PORTAL - PREREQUISITE CHECK IN interface. The top navigation bar includes the portal logo, home icon, and links for My Health Record, My Blood Donations, and Health Centre. A red button labeled "VC APPOINTMENT" is visible in the top right corner. The main content area is titled "Appointment" and "Booking". A progress bar at the top indicates the current step: "Self-monitoring" (step 3) is active, with "Consent" (step 1) and "Patient Instructions" (step 2) completed, and "Done" (step 4) pending.

The "Self-monitoring" section contains a warning: "Self-monitoring is not mandatory to proceed with the virtual consultation session". Below this, there are two dropdown menus for "Vital Sign" and "BMI".

The "GLUCOSE" section is expanded, showing a "GLUCOSE PANEL" with a red circle containing the number "1". It includes two input fields for glucose levels:

- Fasting Glucose: Input field contains "3". Normal Range: 3.9 - 7 mmol/L.
- Random Glucose: Input field contains "3". Normal Range: 3.9 - 11 mmol/L.

Below the input fields is a red "SAVE" button with a red circle containing the number "2". At the bottom of the form are three buttons: "Cancel", "Back", and "Next".

The footer of the page shows "RekodSaya Portal © 2023" and a "Session Timeout in 14 Min 59 Sec" warning.

10. Self-monitoring - Final Glucose

1. The system displays pop-up message.
2. Click on [OK] button.

The screenshot displays the REKODSAYA PORTAL interface during a self-monitoring session. The top navigation bar includes the portal logo, home icon, and menu items: Home, My Health Record, My Blood Donations, and Health Centre. A red button labeled 'VC APPOINTMENT' is visible on the right. The user's name 'KASIH BALQIS' is shown in the top right corner. The main content area is divided into 'Appointment' and 'Booking' sections. The 'Appointment' section shows a progress bar with four steps: Consent, Patient Instructions, Self-monitoring (current step), and Done. Below the progress bar, a message states: 'Self-monitoring is not mandatory to proceed with the virtual consultation session'. The 'GLUCOSE' section contains a 'GLUCOSE PANEL' with two input fields: 'Fasting Glucose' (value: 3, Normal Range: 3.9 - 7 mmol/L) and 'Random Glucose' (value: 3, Normal Range: 3.9 - 11 mmol/L). A red 'SAVE' button is located below the input fields. A pop-up message is displayed in the center, titled 'FINAL GLUCOSE' with a warning icon. The message text reads: 'Are you sure you want to SAVE this GLUCOSE? This action is not be able to update once saved'. The pop-up has 'Cancel' and 'OK' buttons. A blue circle with the number '1' is placed over the 'FINAL GLUCOSE' title, and another blue circle with the number '2' is placed over the 'OK' button. The bottom of the screen shows the footer: 'RekodSaya Portal © 2023' and a session timeout warning: 'Session Timeout in 14 Min 59 Sec'.

11. Submit Prerequisites

1. Click on [Next] button to submit prerequisites.

The screenshot displays the RekodSaya Portal interface during a prerequisite check-in. The top navigation bar includes the portal logo, home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. On the right, there are language options (EN), a notification bell with a '3' badge, a user profile icon, and the name 'KASIH BALQIS'.

The main content area is titled 'Appointment' and 'Booking'. A progress bar at the top shows four steps: 'Consent' (completed), 'Patient Instructions' (completed), 'Self-monitoring' (current step, indicated by a blue circle with '3'), and 'Done' (indicated by a blue circle with '4').

Below the progress bar, a blue notification states: 'Self-monitoring is not mandatory to proceed with the virtual consultation session'. The form contains several input fields:

- 'Vital Sign' dropdown menu.
- 'BMI' dropdown menu.
- 'GLUCOSE' section with a 'GLUCOSE PANEL' sub-section.
- 'Fasting Glucose' input field with the value '3' and a normal range of '3.9 - 7 mmol/L'.
- 'Random Glucose' input field with the value '3' and a normal range of '3.9 - 11 mmol/L'.
- A 'SAVE' button with a camera icon.

At the bottom of the form, there are four buttons: 'Cancel' (red), 'Back' (blue), 'Next' (blue), and a large blue button with the number '1'. The footer of the page shows 'RekodSaya Portal © 2023' and a session timeout warning: 'Session Timeout in 14 Min 59 Sec'.

12. Submit Prerequisites (Continued)

1. The system displays pop-up message.
2. Click on [OK] button.

The screenshot displays the REKODSAYA PORTAL interface. At the top, there is a navigation bar with the portal logo, home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. On the right side of the navigation bar, there is a 'VC APPOINTMENT' button, language settings (EN), a notification bell, and the text 'KASIH BALQIS'.

The main content area is titled 'Appointment' and 'Booking'. A warning message reads: 'If you feel unwell, kindly see your health care provider physically for a more thorough exam'. Below this is the 'PREREQUISITES FORM' with a progress bar showing four steps: 'Consent', 'Self-monitoring', and 'Done'. The 'Consent' step is currently active.

A pop-up message titled 'Submit Prerequisites' (marked with a blue '1') is displayed in the center. The message asks: 'Are you sure you want to SUBMIT this prerequisites? This action is not be able to update once saved'. The pop-up has 'Cancel' and 'OK' buttons. The 'OK' button is highlighted with a blue '2', indicating the next step in the process.

At the bottom of the page, there is a footer with 'RekodSaya Portal © 2023' and a 'Session Timeout in 14 Min 58 Sec' indicator.

13. Info (Self-monitoring)

1. The system displays pop-up message.
2. Click on [OK] button.

The screenshot displays the REKODSAYA PORTAL interface. The top navigation bar includes the portal logo, a home icon, and links for 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre'. The user's name 'KASIH BALQIS' is visible in the top right corner. A left-hand menu lists various health services: Health Events, Appointment, Allergy, Imaging, Family History, Comorbidities, Medication, Laboratory, Observation, Immunization, Transfusion History, Referral, Self Assessment, and Billing. The main content area is titled 'Appointment' and shows a table of bookings. The table is divided into three sections: 'BOOKED', 'FULFILLED', and 'CANCELLED/EXPIRED'. The 'BOOKED' section contains one entry with status 'ARRIVED' and a 'CANCEL' button. The 'FULFILLED' section is currently empty. The 'CANCELLED/EXPIRED' section contains one entry with status 'CANCELLED'. A white pop-up message box is overlaid on the table, containing a warning icon, the text 'INFO', and the message 'You can view self-monitoring entry at observation menu.' The pop-up has an 'OK' button and is annotated with a blue circle containing the number '1' and another blue circle containing the number '2' pointing to the 'OK' button. At the bottom of the page, there is a footer with 'RekodSaya Portal © 2023' and a 'Session Timeout in 14 Min 58 Sec' indicator.

No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	Option
BOOKED									
1	HTAN - KLINIK VC INDIVIDUAL	25-Aug-2023	16:00 - 17:00	BILIK SAKIT 1	Addiction & Substance Abuse	Virtual	Hospital Tuanku Ampuan Najihah, Kuala Pilah	ARRIVED	CANCEL
FULFILLED									
No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	
CANCELLED/EXPIRED									
No.	Schedule Name	Date	Time	Location	Service	Category	Facility	Status	
1	HTAN - KLINIK VC INDIVIDUAL	25-Aug-2023	16:00 - 17:00	BILIK SAKIT 1	Addiction & Substance Abuse	Virtual	Hospital Tuanku Ampuan Najihah, Kuala Pilah	CANCELLED	

14. Self-monitoring - Observation

1. Click on the observation menu button to open observation page.

The screenshot displays the 'Observation' page in the REKODSAYA PORTAL. The page features a navigation menu on the left with 'Observation' selected. The main content area is titled 'Observation' and contains a 'Vital Signs' table. The table has columns for 'No.', 'Activity', 'Normal Range', and a date/time entry '25-Aug-2023 15:40'. The table contains 7 rows of data for various vital signs. The 'Observation' menu item is highlighted in the left sidebar.

No.	Activity	Normal Range	25-Aug-2023 15:40
1	Respiratory Rate	12 - 20 breaths/min	67
2	Heart Rate	60 - 100 beats/min	99
3	Oxygen Saturation	95% - 100%	89
4	Body Temperature	35°C - 37.5°C	37
5	Systolic Blood Pressure	90 - 120 mmHg	66
6	Dystolic Blood Pressure	60 - 80 mmHg	55
7	Pain Score		

15. Self-monitoring - Observation (Select Panel)

1. Choose Vital Signs/ Glucose/ Anthropometry from [Select Panel] dropdown list.

Note:

Vital Signs displays as the default panel

The screenshot displays the 'Observation' section of the REKODSAYA PORTAL. A 'Select Panel' dropdown menu is open, showing the following options:

- Vital Signs
- Glucose
- Anthropometry
- Anthropometry Child
- Dental Chart

The main content area shows a table of vital signs data for the date 25-Aug-2023 15:40. The table has the following structure:

No.	Activity	Normal Range	25-Aug-2023 15:40
1	Respiratory Rate	12 - 20 breaths/min	67
2	Heart Rate	60 - 100 beats/min	99
3	Oxygen Saturation	95% - 100%	89
4	Body Temperature	35°C - 37.5°C	37
5	Systolic Blood Pressure	90 - 120 mmHg	66
6	Dystolic Blood Pressure	60 - 80 mmHg	55
7	Pain Score		

The interface also includes a left sidebar menu with 'Observation' selected, a top navigation bar with 'Home', 'My Health Record', 'My Blood Donations', and 'Health Centre', and a footer with 'RekodSaya Portal © 2023' and 'Session Timeout in 14 Min 57 Sec'.

16. Self-monitoring - Observation (Vital Signs)

1. Mouse hover to [information] icon.

Note:

The system displays details information with recorder name mentioned as 'Patient Self-Recorded' referring to patient self-input.

The screenshot displays the 'Observation' section of the REKODSAYA PORTAL. The page features a navigation menu on the left with 'Observation' selected. The main content area shows a table of Vital Signs recorded on 25-Aug-2023 at 15:40. The table includes columns for No., Activity, Normal Range, and a date/time column. The recorder is identified as 'Patient Self-Recorded'.

No.	Activity	Normal Range	25-Aug-2023 15:40
1	Respiratory Rate	12 - 20 breaths/min	67
2	Heart Rate	60 - 100 beats/min	99
3	Oxygen Saturation	95% - 100%	89
4	Body Temperature	35°C - 37.5°C	37
5	Systolic Blood Pressure	90 - 120 mmHg	66
6	Dystolic Blood Pressure	60 - 80 mmHg	55
7	Pain Score		

At the bottom of the page, there is a footer with 'RekodSaya Portal © 2023' and a session timeout warning: 'Session Timeout in 14 Min 59 Sec'.

17. Self-monitoring - Observation (Glucose)

1. Mouse hover to [information] icon.

Note:

The system displays details information with recorder name mentioned as 'Patient Self-Recorded' referring to patient self-input.

The screenshot displays the RekodSaya Portal interface. At the top, there is a navigation bar with icons for Home, My Health Record, My Blood Donations, and Health Centre. The user is logged in as KASIH BALQIS. A sidebar menu on the left lists various health-related categories, with 'Observation' selected. The main content area shows a table of vital signs and two detailed observation sections: 'Glucose' and 'Anthropometry'.

No.	Activity	Normal Range	Value
1	Respiratory Rate	12 - 20 breaths/min	67
2	Heart Rate	60 - 100 beats/min	99
3	Oxygen Saturation	95% - 100%	89
4	Body Temperature	35°C - 37.5°C	37
5	Systolic Blood Pressure	90 - 120 mmHg	66
6	Dystolic Blood Pressure	60 - 80 mmHg	55
7	Pain Score		

No.	Activity	Normal Range	Value
1	Fasting Glucose	3.9 - 7 mmol/L	3
2	Random Glucose	3.9 - 11 mmol/L	3

The 'Glucose' section includes a date and time stamp '25-Aug-2023 15:41' and a recorder name 'Patient Self-Recorded'. The 'Anthropometry' section also shows a recorder name 'Patient Self-Recorded'.

At the bottom of the page, there is a footer with 'RekodSaya Portal © 2023' and a session timeout warning: 'Session Timeout in 14 Min 59 Sec'.

18. Self-monitoring - Observation (Anthropometry)

1. Mouse hover to [information] icon.

Note:

The system displays details information with recorder name mentioned as 'Patient Self-Recorded' referring to patient self-input

The screenshot displays the 'Observation' section of the RekodSaya Portal. It features a navigation menu on the left and a main content area with two data tables. The top table, titled 'Glucose', shows two entries: 'Fasting Glucose' and 'Random Glucose', both recorded on 25-Aug-2023 at 15:41 with a value of 3. The bottom table, titled 'Anthropometry', shows four entries: 'Height (cm)' (150), 'Weight (kg)' (55), 'BMI (kg/m²)' (24.44), and 'Waist Circumference (cm)' (42), all recorded on the same date and time. Both tables include a 'Patient Self-Recorded' status indicator and an information icon.

Glucose

No.	Activity	Normal Range	25-Aug-2023 15:41
1	Fasting Glucose	3.9 - 7 mmol/L	3
2	Random Glucose	3.9 - 11 mmol/L	3

Anthropometry

No.	Activity	25-Aug-2023 15:41
1	Height (cm)	150
2	Weight (kg)	55
3	BMI (kg/m ²)	24.44
4	Waist Circumference (cm)	42

RekodSaya Portal © 2023 Session Timeout in 14 Min 59 Sec